

MINIMUM SEEDS REQUIRED FOR TESTING



The following are the minimum number of seeds reported for each test. One full sample bag is requested for most testing which allows enough seed for sample retention and retesting, unless otherwise indicated.

TEST	NUMBER OF SEEDS
Warm Germination	400
Sand Germination	400
Cold Test	200
Saturated Cold Test	200
TZ Test	100
Accelerated Aging	200
BioAssay Corn	400
BioAssay Soybeans	200
Elisa (Trait) Test	90
Lateral Flow Strip Test	100—3,000
PCR*	3,000
IsoElectric Focusing	90-and 45 seeds of each parent

*If requesting PCR in addition to other testing, please send 2 sample bags banded together